




performance^x
LINE **PRO-
DUCT
CATA-
LOG**

TRI® Performance Line


Implant System



Narrow, Bone- and Tissue-Level

The TRI® Performance Line is incorporated in the Bone-Level (TRI®-Narrow, TRI®-Vent) and Tissue-Level (TRI®-Octa) implant lines and ensures optimal esthetic results both in the anterior and posterior area.

Lean & intelligent implant concept.



The TRI® Performance Line combines the pure simplicity of one connection with the esthetic flexibility of the TRI® Soft Tissue Concept. A unique lean, and yet intelligent blend for success in daily practice.



Surgical Kit

The TRI® Performance Line is characterized by one unique surgical kit for bone- and tissue-level line. Its intuitive and elegant design prioritizes efficiency and ease of use.

TRI®-Vent Implants - Includes Surgical Cover Screw TVSCS

| Implant Ø | Platform | Connection | 6.5mmL | 8mmL | 10mmL | 11.5mmL | 13mmL | 16mmL |
|-----------|----------|--------------------------|---------|---------|---------|---------|---------|---------|
| 3.75 mm | 3.5 mmD | 2.5 mmD Internal Hexagon | | TV37Mo8 | TV37M10 | TV37M11 | TV37M13 | TV37M16 |
| 4.1 mm | | | TV41Mo6 | TV41Mo8 | TV41M10 | TV41M11 | TV41M13 | TV41M16 |
| 4.7 mm | | | TV47Mo6 | TV47Mo8 | TV47M10 | TV47M11 | TV47M13 | TV47M16 |

TRI®-Narrow Implants - Includes Surgical Cover Screw TNSCS

| Implant Ø | Platform | Connection | 10mmL | 11mmL | 13mmL | 16mmL |
|-----------|----------|--------------------------|---------|---------|---------|---------|
| 3.3 mm | 3.0 mmD | 2.1 mmD Internal Hexagon | TN33M10 | TN33M11 | TN33M13 | TN33M16 |

TRI®-Octa Implants - Includes Surgical Cover Screw TOSCS



| Implant Ø | Platform | Connection | 6.5mmL | 8mmL | 10mmL | 11.5mmL | 13mmL |
|-----------|------------------------|--------------------------|---------|---------|---------|---------|---------|
| 3.75 mm | 4.8 mmD 1.8mm Shoulder | 3.1 mmD Internal Octagon | | TO37Mo8 | TO37M10 | TO37M11 | TO37M13 |
| 4.1 mm | | | TO41Mo6 | TO41Mo8 | TO41M10 | TO41M11 | TO41M13 |
| 4.7 mm | | | TO47Mo6 | TO47Mo8 | TO47M10 | TO47M11 | TO47M13 |

TRI®-Vent Prosthetics


3.5 mmD Implant Platform

HEALING COLLARS

Tissue healing - Healing Collars, sterile



| | Product Code | | Ø | GH | Material | Qty |
|---|--------------|----------------|--------|--------|-----------|-----|
|  | TVHC35-30 | Healing Collar | 3.5 mm | 3.0 mm | Ti-6Al-4V | 1 |
|  | TVHC35-45 | Healing Collar | 3.5 mm | 4.5 mm | Ti-6Al-4V | 1 |
| | TVHC35-60 | Healing Collar | 3.5 mm | 6.0 mm | Ti-6Al-4V | 1 |
| | TVHC50-30 | Healing Collar | 5.0 mm | 3.0 mm | Ti-6Al-4V | 1 |
| | TVHC50-45 | Healing Collar | 5.0 mm | 4.5 mm | Ti-6Al-4V | 1 |
| | TVHC50-60 | Healing Collar | 5.0 mm | 6.0 mm | Ti-6Al-4V | 1 |

Surgical cover screws - sterile

| | Product Code | | Ø | GH | Material | Qty |
|---|--------------|----------------------------|--------|----|-----------|-----|
|  | TVSCS | Surgical Cover Screw | 3.5 mm | | Ti-6Al-4V | 1 |
| | TVSCS-SL | Cover Screw for Sinus Lift | 5.2 mm | | Ti-6Al-4V | 1 |


IMPRESSION TRANSFER COMPONENTS, IMPLANT ANALOGS

Direct transfer for open-tray impressions TRI® 3D-Touch Impression - includes screw DTRS-TV05

| | Product Code | | Ø | GH | Material | Qty |
|---|--------------|---------------------------------|--------|----|-----------|-----|
|  | TV05-35 | Direct Transfer Component | 3.5 mm | | Ti-6Al-4V | 1 |
| | TV05-50 | Direct Transfer Component | 5.0 mm | | Ti-6Al-4V | 1 |
|  | DTRS-TV05 | Direct Transfer Retaining Screw | | | Ti-6Al-4V | 1 |

Indirect transfer/straight abutment for closed-tray impressions

TRI® 2in1 Impression Abutments - includes screws RS-TV10 and RS-TV10-ITC

| | Product Code | | Ø | GH | Material | Qty |
|---|--------------|---|--------|----|-----------|-----|
|  | TV10-35-ITC | Indirect Transfer Component | 3.5 mm | | Ti-6Al-4V | 1 |
| | TV10-50-ITC | Indirect Transfer Component | 5.0 mm | | Ti-6Al-4V | 1 |
| | RS-TV10-ITC | Indirect Transfer Component Impression Screw | | | POM | 1 |
| | RS-TV10 | Replacement Retaining Screw | | | Ti-6Al-4V | 1 |
| | RS-TV10-Lab | Replacement Retaining Screw - Lab Use - Green | | | Ti-6Al-4V | 1 |

GH = Gingiva Height

Implant Analog



| Product Code | | Material | Qty |
|--------------|---------------------------------|----------|-----|
| TVA-1 | Digital Implant Analog | Steel | 1 |
| TVA-5 | Digital Implant Analog (5 Pack) | Steel | 5 |

COMPONENTS FOR CEMENT-RETAINED RESTORATIONS

Plastic Temporary Abutments - includes screw RS-TV10/RS-TV20



| Product Code | | Ø | Material | Qty |
|--------------|---|--------|-----------|-----|
| TV10-PTA | Temporary Abutment, Straight | 4.5 mm | PEEK | 1 |
| RS-TV20 | Replacement Retaining Screw | | Ti-6Al-4V | 1 |
| RS-TV20-Lab | Replacement Retaining Screw - Lab Use - Red | | Ti-6Al-4V | 1 |
| TV10-EHA-50* | Esthetic Healing Abutment | 5.0 mm | PEEK | 1 |
| TV10-EHA-70* | Esthetic Healing Abutment | 7.0 mm | PEEK | 1 |
| RS-TV10 | Replacement Retaining Screw | | Ti-6Al-4V | 1 |
| RS-TV10-Lab | Replacement Retaining Screw - Lab Use - Green | | Ti-6Al-4V | 1 |

Straight Abutments with TRI®-Friction - includes screw RS-TV10



| Product Code | | Ø | GH | Material | Qty |
|--------------|---|--------|--------|-----------|-----|
| TV10-35-F | Straight Abutment with TRI®-Friction | 3.5 mm | | Ti-6Al-4V | 1 |
| TV10-50-F | Straight Abutment with TRI®-Friction | 5.0 mm | | Ti-6Al-4V | 1 |
| TV10-41-F | Contour Abutment, Straight with TRI®-Friction | 4.5 mm | 1.0 mm | Ti-6Al-4V | 1 |
| TV10-43-F | Contour Abutment, Straight with TRI®-Friction | 4.5 mm | 3.0 mm | Ti-6Al-4V | 1 |
| RS-TV10 | Replacement Retaining Screw | | | Ti-6Al-4V | 1 |
| RS-TV10-Lab | Replacement Retaining Screw - Lab Use - Green | | | Ti-6Al-4V | 1 |

Angled Abutments with TRI®-Friction - includes screws RS-TV10/RS-TV20




| Product Code | | Ø | GH | Material | Qty |
|--------------|--|--------|--------|-----------|-----|
| TV20-50-F | Angled Abutment, Standard - 20° with TRI®-Friction | 5.0 mm | | Ti-6Al-4V | 1 |
| RS-TV20 | Replacement Retaining Screw | | | Ti-6Al-4V | 1 |
| RS-TV20-Lab | Replacement Retaining Screw - Lab Use - Red | | | Ti-6Al-4V | 1 |
| TV20-41-F | Contour Abutment, Angled - 17° with TRI®-Friction | 4.5 mm | 1.0 mm | Ti-6Al-4V | 1 |
| TV20-43-F | Contour Abutment, Angled - 17° with TRI®-Friction | 4.5 mm | 3.0 mm | Ti-6Al-4V | 1 |
| RS-TV10 | Replacement Retaining Screw | | | Ti-6Al-4V | 1 |
| RS-TV10-Lab | Replacement Retaining Screw - Lab Use - Green | | | Ti-6Al-4V | 1 |

*sterile

CASTABLE ABUTMENTS

Castable Abutments, Engaging - includes screw RS-TV10/RS-TV20



| Product Code | | Ø | GH | Material | Qty |
|--------------|---|--------|----|-----------|-----|
| TV30-50-F | Gold Castable Abutment with TRI®-Friction | 5.0 mm | | HMA | 1 |
| RS-TV10 | Replacement Retaining Screw | | | Ti-6Al-4V | 1 |
| RS-TV10-Lab | Replacement Retaining Screw - Lab Use - Green | | | Ti-6Al-4V | 1 |
| TV35-50 | Platin-Iridium Castable Abutment | 5.0 mm | | PtIr | 1 |
| RS-TV20 | Replacement Retaining Screw | | | Ti-6Al-4V | 1 |
| RS-TV20-Lab | Replacement Retaining Screw - Lab Use - Red | | | Ti-6Al-4V | 1 |

COMPONENTS FOR SCREW-RETAINED RESTORATIONS

Screw-retained Multi-Unit Abutments



| Product Code | | Ø | GH | Material | Qty |
|--------------|--|--------|--------|-----------|-----|
| TV40-01* | Screw-Retained Abutment, Straight | 4.5 mm | 1.0 mm | Ti-6Al-4V | 1 |
| TV40-02* | Screw-Retained Abutment, Straight | 4.5 mm | 2.0 mm | Ti-6Al-4V | 1 |
| TV40-04* | Screw-Retained Abutment, Straight | 4.5 mm | 4.0 mm | Ti-6Al-4V | 1 |
| TV40-06* | Screw-Retained Abutment, Straight | 4.5 mm | 6.0 mm | Ti-6Al-4V | 1 |
| TV50-17* | Screw Retained Abutment, 17° Angled | 4.5 mm | 2.0 mm | Ti-6Al-4V | 1 |
| TV50-30* | Screw Retained Abutment, 30° Angled | 4.5 mm | 3.0 mm | Ti-6Al-4V | 1 |
| HC-SRAS* | Healing Cap for Screw-Retained Abutment | | | Ti-6Al-4V | 1 |
| TV40-SCAN | Multi-Unit Scanbody | | | PEEK | 1 |
| TTA | Titanium Temporary Abutment incl. Screw RS-TTA | | | Ti-6Al-4V | 1 |
| RS-TTA | Replacement Retaining Screw for TTA | | | Ti-6Al-4V | 1 |
| RS-TTAL | Replacement Retaining Screw Long for TTA | | | Ti-6Al-4V | 1 |
| DTC-SRA | Direct Transfer Component for Screw-Retained Abutment incl. DTRS-SRA | | | Ti-6Al-4V | 1 |
| DTRS-SRA | Direct Transfer Retaining Screw | | | Ti-6Al-4V | 1 |
| MUA | Multi-Unit Digital Implant Analog | | | Ti-6Al-4V | 1 |
| MUA-5 | Multi-Unit Digital Implant Analog | | | Ti-6Al-4V | 5 |
| RS-MUA-D-4 | Replacement Retaining Screw for MUA | | | Ti-6Al-4V | 4 |
| PCC | Plastic Castable Coping incl. Screw RS-PCC | | | POM | 1 |
| RS-PCC | Replacement Retaining Screw for PCC | | | Ti-6Al-4V | 1 |
| RS-PCC-Lab | Replacement Retaining Screw - Lab Use - Blue | | | Ti-6Al-4V | 1 |

*sterile

OVERDENTURE RESTORATIONS

Ball Abutments - includes CAHL



| Product Code | | Ø | GH | Material | Qty |
|--------------|---|--------|--------|-----------------|-----|
| TV60-02 | Ball Abutment incl. CAHL | 4.5 mm | 2.0 mm | Ti-6Al-4V | 1 |
| TV60-04 | Ball Abutment incl. CAHL | 4.5 mm | 4.0 mm | Ti-6Al-4V | 1 |
| CAHL | Cap Attachment Housing and Liner - Transparent | | | Stainless steel | 1 |
| CAH | Cap Attachment Housing | | | Stainless steel | 1 |
| CALT | Cap Attachment Liner - (Transparent - Standard retention) | | | Nylon | 1 |
| CALG | Cap Attachment Liner - (Gray - Rigid retention) | | | Nylon | 1 |
| CALB | Cap Attachment Liner - Lab Use - Black | | | Nylon | 1 |


TRI® Raptor Abutments - includes Caps Assortment Kit (Prod. Code 192ECE)



| Product Code | | Ø | GH | Material | Qty |
|--------------|---|---|--------|-----------|-----|
| 130CV.01 | TRI® Raptor Abutment | | 1.0 mm | TiN | 1 |
| 130CV.02 | TRI® Raptor Abutment | | 2.0 mm | TiN | 1 |
| 130CV.03 | TRI® Raptor Abutment | | 3.0 mm | TiN | 1 |
| 130CV.04 | TRI® Raptor Abutment | | 4.0 mm | TiN | 1 |
| 130CV.05 | TRI® Raptor Abutment | | 5.0 mm | TiN | 1 |
| 130CV.06 | TRI® Raptor Abutment | | 6.0 mm | TiN | 1 |
| 140CEV | Retentive Caps/violet (strong retention 2700g) | | | Nylon | 4 |
| 140CET | Retentive Caps/clear (standard retention 1800g) | | | Nylon | 4 |
| 140CER | Retentive Caps/pink (soft retention 1200g) | | | Nylon | 4 |
| 140CEG | Retentive Caps/yellow (extra-soft retention 600g) | | | Nylon | 4 |
| 192ECE | Caps Assortment Kit | | | Ti, Nylon | 1 |
| 335SBC | Smart Box Kit | | | Ti, Nylon | 1 |
| 144AE | Analog for Lab | | | Stainless | 2 |
| 044CAIN | Transfer for Impression | | | Stainless | 2 |

COMPONENTS FOR TRI+ DIGITAL CAD/CAM

CAD/CAM Abutments - includes screw RS-TV10



| Product Code | | Ø | GH | Material | Qty |
|--------------|---|--------|--------|-----------|-----|
| TV70-SCAN | Scanbody for Laboratory and Intraoral Scanner | | | PEEK | 1 |
| TV70-07-F | Engaging Titanium Bonding-Base with TRI®-Friction | 4.5 mm | 0.7 mm | Ti-6Al-4V | 1 |
| TV70-20-F | Engaging Titanium Bonding-Base with TRI®-Friction | 4.5 mm | 2.0 mm | Ti-6Al-4V | 1 |
| TV70-07-NE | Non-Engaging Titanium Bonding-Base | 4.5 mm | 0.7 mm | Ti-6Al-4V | 1 |
| TV70-20-NE | Non-Engaging Titanium Bonding-Base | 4.5 mm | 2.0 mm | Ti-6Al-4V | 1 |
| TV70-Blank-F | Titanium Milling Blank with TRI®-Friction | | | Ti-6Al-4V | 1 |
| RS-TV10 | Replacement Retaining Screw | | | Ti-6Al-4V | 1 |
| RS-TV10-Lab | Replacement Retaining Screw - Lab Use - Green | | | Ti-6Al-4V | 1 |

GH = Gingiva Height

TRI®-Narrow Prosthetics

3.0 mmD Implant Platform

HEALING COLLARS

Tissue healing - Healing Collars, sterile



TNHC35-30

| Product Code | | Ø | GH | Material | Qty |
|--------------|----------------|--------|--------|-----------|-----|
| TNHC35-30 | Healing Collar | 3.5 mm | 3.0 mm | Ti-6Al-4V | 1 |
| TNHC35-45 | Healing Collar | 3.5 mm | 4.5 mm | Ti-6Al-4V | 1 |
| TNHC35-60 | Healing Collar | 3.5 mm | 6.0 mm | Ti-6Al-4V | 1 |
| TNHC45-30 | Healing Collar | 4.5 mm | 3.0 mm | Ti-6Al-4V | 1 |
| TNHC45-50 | Healing Collar | 4.5 mm | 5.0 mm | Ti-6Al-4V | 1 |

Surgical cover screws, sterile



TNSCS

| Product Code | | Ø | GH | Material | Qty |
|--------------|----------------------|--------|----|-----------|-----|
| TNSCS | Surgical Cover Screw | 3.2 mm | | Ti-6Al-4V | 1 |

IMPRESSION TRANSFER COMPONENTS, IMPLANT ANALOGS

Direct transfer for open-tray impressions - TRI® 3D-Touch Impression - includes screw DTRS-TN05



TN05-35



DTRS-TN05

| Product Code | | Ø | GH | Material | Qty |
|--------------|---------------------------------|--------|----|-----------|-----|
| TN05-35 | Direct Transfer Component | 3.5 mm | | Ti-6Al-4V | 1 |
| DTRS-TN05 | Direct Transfer Retaining Screw | | | Ti-6Al-4V | 1 |

Indirect transfer/straight abutment for closed-tray impressions

TRI® 2in1 Impression Abutments - includes screw RS-TN10



TN 10-35-ITC

| Product Code | | Ø | GH | Material | Qty |
|--------------|-----------------------------|--------|----|-----------|-----|
| TN10-35-ITC | Indirect Transfer Component | 3.5 mm | | Ti-6Al-4V | 1 |
| RS-TN10 | Replacement Retaining Screw | | | Ti-6Al-4V | 1 |

Implant Analog



TNA



TNIA

| Product Code | | Material | Qty |
|--------------|---------------------------------|-----------|-----|
| TNA-1 | Digital Implant Analog | Ti-6Al-4V | 1 |
| TNA-5 | Digital Implant Analog (5 Pack) | Ti-6Al-4V | 5 |

GH = Gingiva Height

COMPONENTS FOR CEMENT-RETAINED RESTORATIONS

Plastic Temporary Abutments - includes screw RS-TN20



TN10-PTA

| Product Code | | Ø | Material | Qty |
|--------------|---|--------|-----------|-----|
| TN10-PTA | Temporary Abutment, Straight | 3.5 mm | PEEK | 1 |
| RS-TN20 | Replacement Retaining Screw | | Ti-6Al-4V | 1 |
| RS-TN20-Lab | Replacement Retaining Screw - Lab Use - Red | | Ti-6Al-4V | 1 |

Straight Abutments with TRI®-Friction - includes screw RS-TN10



TN10-41-F



TN10-35-F

| Product Code | | Ø | GH | Material | Qty |
|--------------|---|--------|--------|-----------|-----|
| TN10-35-F | Straight Abutment with TRI®-Friction | 3.5 mm | | Ti-6Al-4V | 1 |
| TN10-41-F | Contour Abutment, Straight with TRI®-Friction | 4.0 mm | 1.0 mm | Ti-6Al-4V | 1 |
| TN10-43-F | Contour Abutment, Straight with TRI®-Friction | 4.0 mm | 3.0 mm | Ti-6Al-4V | 1 |
| RS-TN10 | Replacement Retaining Screw | | | Ti-6Al-4V | 1 |
| RS-TN10-Lab | Replacement Retaining Screw - Lab Use - Green | | | Ti-6Al-4V | 1 |

Angled Abutments with TRI®-Friction - includes screws RS-TN10



TN20-41-F



TN20-45-F

| Product Code | | Ø | GH | Material | Qty |
|--------------|--|--------|--------|-----------|-----|
| TN20-45-F | Angled Abutment, Standard - 20° with TRI®-Friction | 4.5 mm | | Ti-6Al-4V | 1 |
| TN20-41-F | Contour Abutment, Angled - 17° with TRI®-Friction | 4.0 mm | 1.0 mm | Ti-6Al-4V | 1 |
| TN20-43-F | Contour Abutment, Angled - 17° with TRI®-Friction | 4.0 mm | 3.0 mm | Ti-6Al-4V | 1 |
| RS-TN10 | Replacement Retaining Screw | | | Ti-6Al-4V | 1 |
| RS-TN10-Lab | Replacement Retaining Screw - Lab Use - Green | | | Ti-6Al-4V | 1 |

CASTABLE ABUTMENTS

Gold abutments, engaging - includes screw RS-TN10



TN30-45-F

| Product Code | | Ø | Material | Qty |
|--------------|---|--------|-----------|-----|
| TN30-45-F | Gold Castable Abutment with TRI®-Friction | 4.5 mm | HMA | 1 |
| RS-TN10 | Replacement Retaining Screw | | Ti-6Al-4V | 1 |
| RS-TN10-Lab | Replacement Retaining Screw - Lab Use - Green | | Ti-6Al-4V | 1 |

COMPONENTS FOR TRI+ DIGITAL CAD/CAM

CAD/CAM Abutments - includes screw RS-TN10



TN70-SCAN



TN70-07-F



TN70-07-NE



TN70-Blank-F

| Product Code | | Ø | GH | Material | Qty |
|--------------|---|--------|--------|-----------|-----|
| TN70-SCAN | Scanbody for Laboratory and Intraoral Scanner | | | PEEK | 1 |
| TN70-07-F | Engaging Titanium Bonding-Base with TRI®-Friction | 4.2 mm | 0.7 mm | Ti-6Al-4V | 1 |
| TN70-07-NE | Non-Engaging Titanium Bonding-Base | 4.2 mm | 0.7 mm | Ti-6Al-4V | 1 |
| TN70-Blank-F | Titanium Milling Blank with TRI®-Friction | | | Ti-6Al-4V | 1 |
| RS-TN10 | Replacement Retaining Screw | | | Ti-6Al-4V | 1 |
| RS-TN10-Lab | Replacement Retaining Screw - Lab Use - Green | | | Ti-6Al-4V | 1 |

COMPONENTS FOR SCREW-RETAINED RESTORATIONS

Screw-retained Multi-Unit Abutments



| Product Code | | Ø | GH | Material | Qty |
|--------------|--|--------|--------|-----------|-----|
| TN40-02* | Screw-Retained Abutment, Straight | 4.5 mm | 2.0 mm | Ti-6Al-4V | 1 |
| TN40-04* | Screw-Retained Abutment, Straight | 4.5 mm | 4.0 mm | Ti-6Al-4V | 1 |
| HC-SRAS* | Healing Cap for Screw-Retained Abutment | | | Ti-6Al-4V | 1 |
| TV40-SCAN | Multi-Unit Scanbody | | | PEEK | 1 |
| TTA | Titanium Abutment incl. Screw RS-TTA | | | Ti-6Al-4V | 1 |
| RS-TTA | Replacement Retaining Screw for TTA | | | Ti-6Al-4V | 1 |
| RS-TTAL | Replacement Retaining Screw Long for TTA | | | Ti-6Al-4V | 1 |
| DTC-SRA | Direct Transfer Component for Screw-Retained Abutment incl. DTRS-SRA | | | Ti-6Al-4V | 1 |
| DTRS-SRA | Direct Transfer Retaining Screw | | | Ti-6Al-4V | 1 |
| MUA | Multi-Unit Digital Implant Analog | | | Ti-6Al-4V | 1 |
| MUA-5 | Multi-Unit Digital Implant Analog | | | Ti-6Al-4V | 5 |
| RS-MUA-D-4 | Replacement Retaining Screw for MUA | | | Ti-6Al-4V | 4 |
| PCC | Plastic Castable Coping incl. Screw RS-PCC | | | POM | 1 |
| RS-PCC | Replacement Retaining Screw for PCC | | | Ti-6Al-4V | 1 |
| RS-PCC-Lab | Replacement Retaining Screw - Lab Use - Blue | | | Ti-6Al-4V | 1 |

OVERDENTURE RESTORATIONS

Ball Abutments - includes CAHL



| Product Code | | Ø | GH | Material | Qty |
|--------------|---|--------|--------|-----------------|-----|
| TN60-02 | Ball Abutment incl. CAHL | 4.5 mm | 2.0 mm | Ti-6Al-4V | 1 |
| TN60-04 | Ball Abutment incl. CAHL | 4.5 mm | 4.0 mm | Ti-6Al-4V | 1 |
| CAHL | Cap Attachment Housing and Liner - Transparent | | | Stainless steel | 1 |
| CAH | Cap Attachment Housing | | | Stainless steel | 1 |
| CALT | Cap Attachment Liner - (Transparent - Standard retention) | | | Nylon | 1 |
| CALG | Cap Attachment Liner - (Gray - Rigid retention) | | | Nylon | 1 |
| CALB | Cap Attachment Liner - Lab Use - Black | | | Nylon | 1 |

TRI® Raptor Abutments - includes Caps Assortment Kit (Prod. Code 192ECE)



| Product Code | | Ø | GH | Material | Qty |
|--------------|---|---|--------|-----------|-----|
| 130TN.01 | TRI® Raptor Abutment | | 1.0 mm | TiN | 1 |
| 130TN.02 | TRI® Raptor Abutment | | 2.0 mm | TiN | 1 |
| 130TN.03 | TRI® Raptor Abutment | | 3.0 mm | TiN | 1 |
| 130TN.04 | TRI® Raptor Abutment | | 4.0 mm | TiN | 1 |
| 130TN.05 | TRI® Raptor Abutment | | 5.0 mm | TiN | 1 |
| 140CEV | Retentive Caps/violet (strong retention 2700g) | | | Nylon | 4 |
| 140CET | Retentive Caps/clear (standard retention 1800g) | | | Nylon | 4 |
| 140CER | Retentive Caps/pink (soft retention 1200g) | | | Nylon | 4 |
| 140CEG | Retentive Caps/yellow (extra-soft retention 600g) | | | Nylon | 4 |
| 192ECE | Caps Assortment Kit | | | Ti, Nylon | 1 |
| 335SBC | Smart Box Kit | | | Ti, Nylon | 1 |
| 144AE | Analog for Lab | | | Stainless | 2 |
| 044CAIN | Transfer for Impression | | | Stainless | 2 |

TRI®-Octa Prosthetics 4.8 mmD Implant Platform

HEALING COLLARS


Tissue healing - Healing Collars, sterile



| Product Code | | Ø | GH | Material | Qty |
|--------------|----------------------|--------|--------|-----------|-----|
| TOSCS | Surgical Cover Screw | 3.5 mm | | Ti-6Al-4V | 1 |
| TOHC51-15 | Healing Collar | 5.1 mm | 1.5 mm | Ti-6Al-4V | 1 |
| TOHC51-30 | Healing Collar | 5.1 mm | 3.0 mm | Ti-6Al-4V | 1 |

IMPRESSION TRANSFER COMPONENTS, IMPLANT ANALOGS

Direct transfer for open-tray impressions - TRI® 3D-Touch Impression - includes screw DTRS-TO05



| Product Code | | Ø | GH | Material | Qty |
|--------------|---------------------------------|--------|----|-----------|-----|
| TO05-51 | Direct Transfer Component | 5.1 mm | | Ti-6Al-4V | 1 |
| DTRS-TO05 | Direct Transfer Retaining Screw | | | Ti-6Al-4V | 1 |

Indirect transfer/straight abutment for closed-tray impressions

TRI® 2in1 Impression Abutments - includes screw RS-TO



TO10-51-ITC

| Product Code | | Ø | GH | Material | Qty |
|--------------|-----------------------------|--------|----|-----------|-----|
| TO10-51-ITC | Indirect Transfer Component | 5.1 mm | | Ti-6Al-4V | 1 |
| RS-TO | Replacement Retaining Screw | | | Ti-6Al-4V | 1 |

Implant Analog



| Product Code | | Material | Qty |
|--------------|---------------------------------|-----------|-----|
| TOA-1 | Digital Implant Analog | Ti-6Al-4V | 1 |
| TOA-5 | Digital Implant Analog (5 Pack) | Ti-6Al-4V | 5 |

GH = Gingiva Height

COMPONENTS FOR CEMENT-RETAINED RESTORATIONS

Plastic Temporary Abutments - includes screw RS-TO



TO10-PTA

| Product Code | | Ø | Material | Qty |
|--------------|------------------------------|--------|-----------|-----|
| TO10-PTA | Temporary Abutment, Straight | 5.1 mm | PEEK | 1 |
| RS-TO | Replacement Retaining Screw | | Ti-6Al-4V | 1 |

Straight Abutments with TRI® Friction - includes screw RS-TO



TO10-35-F



TO10-51-F

| Product Code | | Ø | GH | Material | Qty |
|--------------|--------------------------------------|--------|----|-----------|-----|
| TO10-35-F | Straight Abutment with TRI®-Friction | 3.5 mm | | Ti-6Al-4V | 1 |
| TO10-51-F | Straight Abutment with TRI®-Friction | 5.1 mm | | Ti-6Al-4V | 1 |
| RS-TO | Replacement Retaining Screw | | | Ti-6Al-4V | 1 |

Angled Abutments with TRI®-Friction - includes screws RS-TO



TO20-51-F

| Product Code | | Ø | GH | Material | Qty |
|--------------|--|--------|----|-----------|-----|
| TO20-51-F | Angled Abutment, Standard - 20° with TRI®-Friction | 5.1 mm | | Ti-6Al-4V | 1 |
| RS-TO | Replacement Retaining Screw | | | Ti-6Al-4V | 1 |

COMPONENTS FOR SCREW-RETAINED RESTORATIONS

Screw-retained Multi-Unit Abutments



TO40-16



TO40-30



TTA



MUA



PCC



DTC-SRA



DTRS-SRA

| Product Code | | Ø | GH | Material | Qty |
|--------------|--|--------|--------|-----------|-----|
| TO40-16* | Screw-Retained Abutment, Straight | 4.5 mm | 1.6 mm | Ti-6Al-4V | 1 |
| TO40-30* | Screw-Retained Abutment, Straight | 4.5 mm | 3.0 mm | Ti-6Al-4V | 1 |
| HC-SRAS* | Healing Cap for Screw-Retained Abutment | | | Ti-6Al-4V | 1 |
| TV40-SCAN | Multi-Unit Scanbody | | | PEEK | 1 |
| TTA | Titanium Abutment incl. Screw RS-TTA | | | Ti-6Al-4V | 1 |
| RS-TTA | Replacement Retaining Screw for TTA | | | Ti-6Al-4V | 1 |
| RS-TTAL | Replacement Retaining Screw Long for TTA | | | Ti-6Al-4V | 1 |
| DTC-SRA | Direct Transfer Component for Screw-Retained Abutment incl. DTRS-SRA | | | Ti-6Al-4V | 1 |
| DTRS-SRA | Direct Transfer Retaining Screw | | | Ti-6Al-4V | 1 |
| MUA | Multi-Unit Digital Implant Analog | | | Ti-6Al-4V | 1 |
| MUA-5 | Multi-Unit Digital Implant Analog | | | Ti-6Al-4V | 5 |
| RS-MUA-D-4 | Replacement Retaining Screw for MUA | | | Ti-6Al-4V | 4 |
| PCC | Plastic Castable Coping incl. Screw RS-PCC | | | POM | 1 |
| RS-PCC | Replacement Retaining Screw for PCC | | | Ti-6Al-4V | 1 |
| RS-PCC-Lab | Replacement Retaining Screw - Lab Use - Blue | | | Ti-6Al-4V | 1 |

OVERDENTURE RESTORATIONS

Ball Abutments - includes CAHL



TO60-16



TO60-30



CAHL

| Product Code | | Ø | GH | Material | Qty |
|--------------|---|--------|--------|-----------------|-----|
| To60-16 | Ball Abutment incl. CAHL | 4.5 mm | 1.6 mm | Ti-6Al-4V | 1 |
| TO60-30 | Ball Abutment incl. CAHL | 4.5 mm | 3.0 mm | Ti-6Al-4V | 1 |
| CAHL | Cap Attachment Housing and Liner - Transparent | | | Stainless steel | 1 |
| CAH | Cap Attachment Housing | | | Stainless steel | 1 |
| CALT | Cap Attachment Liner - (Transparent - Standard retention) | | | Nylon | 1 |
| CALG | Cap Attachment Liner - (Gray - Rigid retention) | | | Nylon | 1 |
| CALB | Cap Attachment Liner - Lab Use - Black | | | Nylon | 1 |

TRI® Raptor Abutments - includes Caps Assortment Kit (Prod. Code 192ECE)



144AE



192ECE



044CAIN

| Product Code | | Ø | GH | Material | Qty |
|--------------|---|---|--------|-----------|-----|
| 130BFT.00 | TRI® Raptor Abutment | | 0.0 mm | TiN | 1 |
| 130BFT.01 | TRI® Raptor Abutment | | 1.0 mm | TiN | 1 |
| 130BFT.02 | TRI® Raptor Abutment | | 2.0 mm | TiN | 1 |
| 130BFT.03 | TRI® Raptor Abutment | | 3.0 mm | TiN | 1 |
| 130BFT.04 | TRI® Raptor Abutment | | 4.0 mm | TiN | 1 |
| 130BFT.05 | TRI® Raptor Abutment | | 5.0 mm | TiN | 1 |
| 140CEV | Retentive Caps/violet (strong retention 2700g) | | | Nylon | 4 |
| 140CET | Retentive Caps/clear (standard retention 1800g) | | | Nylon | 4 |
| 140CER | Retentive Caps/pink (soft retention 1200g) | | | Nylon | 4 |
| 140CEG | Retentive Caps/yellow (extra-soft retention 600g) | | | Nylon | 4 |
| 192ECE | Caps Assortment Kit | | | Ti, Nylon | 1 |
| 335SBC | Smart Box Kit | | | Ti, Nylon | 1 |
| 144AE | Analog for Lab | | | Stainless | 2 |
| 044CAIN | Transfer for Impression | | | Stainless | 2 |

COMPONENTS FOR TRI+ DIGITAL CAD/CAM

CAD/CAM Abutments - includes screw RS-TO



TO70-SCAN



TO70-07-F



TO70-07-NE



TO70-Blank-F

| Product Code | | Ø | GH | Material | Qty |
|--------------|---|--------|--------|-----------|-----|
| TO70-SCAN | Scanbody for Laboratory and Intraoral Scanner | | | PEEK | 1 |
| TO70-07-F | Engaging Titanium Bonding-Base with TRI®-Friction | 5.1 mm | 0.7 mm | Ti-6Al-4V | 1 |
| TO70-07-NE | Non-Engaging Titanium Bonding-Base | 5.1 mm | 0.7 mm | Ti-6Al-4V | 1 |
| TO70-Blank-F | Titanium Milling Blank with TRI®-Friction | | | Ti-6Al-4V | 1 |
| RS-TO | Replacement Retaining Screw | | | Ti-6Al-4V | 1 |

GH = Gingiva Height

Starting Drills

| | | |
|--------|-------|----------------|
| LD | Nr. 1 | Lance Drill |
| RB-2.3 | Nr. 2 | Round Bur |
| DE | Nr. 3 | Drill Extender |

Surgical Drills

| | | |
|----------|--------|---|
| TPD-2.3S | Nr. 4 | Pilot Drill - Short - 2.3 mm |
| TPD-2.3L | Nr. 5 | Pilot Drill - Long - 2.3 mm |
| TND-3.3L | Nr. 6 | TRI®-Narrow Final Drill - Long - 2.4/2.8 mm |
| TAP3.3 | Nr. 7 | TRI®-Narrow TAP |
| SBD-3.7S | Nr. 8 | 3.7 Soft Bone Drill - Short - 2.7/3.1 mmD |
| SBD-3.7L | Nr. 9 | 3.7 Soft Bone Drill - Long - 2.7/3.1 mmD |
| DBD-3.7L | Nr. 10 | 3.7 Dense Bone Drill - Long - 2.9/3.6 mmD |
| TAP3.7 | Nr. 11 | 3.7 Tap |
| SBD-4.1S | Nr. 12 | 4.1 Soft Bone Drill - Short - 3.0/3.4 mmD |
| SBD-4.1L | Nr. 13 | 4.1 Soft Bone Drill - Long - 3.0/3.4 mmD |
| DBD-4.1L | Nr. 14 | 4.1 Dense Bone Drill - Long - 3.2/3.9 mmD |
| TAP4.1 | Nr. 15 | 4.1 Tap |
| SBD-4.7S | Nr. 16 | 4.7 Soft Bone Drill - Short - 3.6/4.0 mmD |
| SBD-4.7L | Nr. 17 | 4.7 Soft Bone Drill - Long - 3.6/4.0 mmD |
| DBD-4.7L | Nr. 18 | 4.7 Dense Bone Drill - Long - 3.8/4.6 mmD |
| TAP4.7 | Nr. 19 | 4.7 Tap |
| PP-L | Nr. 21 | Paralell Pin - Long |

Implant Drivers

| | |
|----------|---|
| VD -S | Performance Line TRI®-Vent Implant Driver - Short |
| VD-L | Performance Line TRI®-Vent Implant Driver - Long |
| ND-S | Performance Line TRI®-Narrow Implant Driver - Short |
| ND-L | Performance Line TRI®-Narrow Implant Driver - Long |
| VD-S-TW* | Performance Line TRI®-Vent Implant Driver - Short incl. Torque Wrench |
| VD-L-TW* | Performance Line TRI®-Vent Implant Driver - Long incl. Torque Wrench |
| OD-S | Performance Line TRI®-Octa Implant Driver - Short |
| OD-L | Performance Line TRI®-Octa Implant Driver - Long |
| MD-S-ISO | matrix® Implant Driver - Short |
| MD-L-ISO | matrix® Implant Driver - Long |
| MD-S-TW* | matrix® Implant Driver Short incl. Torque Wrench |
| MD-L-TW* | matrix® Implant Driver Long incl. Torque Wrench |

Torque Wrench

| | |
|--------------|-----------------------|
| TW | Torque Wrench |
| TW - Adapter | Torque Wrench Adapter |

Prosthetic Tools

| | |
|---------------|---|
| PD-L | Prosthetic Driver - Long |
| PD-S | Prosthetic Driver - Short |
| PD-ISO-L* | Prosthetic Driver for Motor - Long |
| PD-ISO-L-Ang* | Prosthetic Driver Angulated - Long |
| TORT | Removal Tool for Performance Line TRI®-Octa Friction Abutments |
| TRT | Removal Tool for Performance Line TRI®-Vent & Narrow Friction Abutments |
| 487.ICE | TRI® Raptor Insertion and Extractor Tool for caps |

Drill Stops


| | |
|--------|------------------------------------|
| TDS6* | Drill Stop / 6 mm |
| TDS8* | Drill Stop / 8 mm |
| TDS10* | Drill Stop / 10 mm |
| TDS11* | Drill Stop / 11 mm |
| TDS13* | Drill Stop / 13 mm |
| TPDS* | Guided Surgery Sleeve for TPD-2.3L |

Other

| | |
|----------|------------------------------|
| TSK | TRI® Surgical Kit - Complete |
| TSK-Tray | TRI® Surgical Kit - Tray |


Prosthetic Kit

| | |
|--------------|---|
| PK | TRI® Prosthetic Kit (includes: PD-ISO-L, PD-L, PD-S, TORT, TRT, TW, TW-Adapter) |
| PK-Tray | Prosthetic Kit-Tray |
| PD-ISO-L* | Prosthetic Driver for Motor - Long |
| PD-L | Prosthetic Driver - Long |
| PD-S | Prosthetic Driver - Short |
| TORT | Removal Tool for Performance Line TRI®-Octa Friction Abutments |
| TRT | Removal Tool for Performance Line TRI®-Vent & Narrow Friction Abutments |
| TW | Torque Wrench |
| TW - Adapter | Torque Wrench Adapter |




* not included in TSK

Drill & Implant Lengths




Drill Sequence


3.3 mmD




3.75 mmD




4.1 mmD



4.7 mmD



Torque Requirements



TRI[®] Pod Handling

3 Options for Implant Insertion



Ratchet



Handpiece



By hand



Step 1
Open Lid



Step 2
Rotate To Remove
implant



Step 3
Remove Cover Screw

TO
TOUCH IT
IS TO
LOVE IT.

The implant packaging by TRI[®] integrates the latest engineering knowhow and innovation. Its revolutionary touchless delivery concept allows to pick-up the implant directly with the surgical hand-piece in one single step. Whilst guaranteeing the highest purity of the TRI[®] SBA surface. A true design innovation which will take your surgical performance to a new level. **Touch it, and you will love it!**

TRIACADEMY

To learn more about our products visit our TRI[®] Academy
Courses and Webinars:
academy@tri-implants.swiss




Become part of
the Community!



Join our
Facebook Group!

DESIGNED AND MANUFACTURED BY TRI[®] IN SWITZERLAND

Designed in the heart of Switzerland, all implants are produced in specialized manufacturing units for dental implants. Therefore benefiting from more than 25 years of medical technology production know-how - the only way to guarantee Swiss Precision with high quality.
Visit us in Switzerland: visit@tri-implants.swiss



TRI Dental Implants Int. AG
Switzerland
00800 3313 3313
www.tri.swiss



TRI® Dental Implants logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of TRI Dental Implants Int. AG.

TRI® Performance Line Catalog EN 2022-12



CE 0297 | FDA Registered